

# CONVERT™ Calf Protocol

CONVERT™ Calf Products are a special blend of naturally occurring direct-fed microorganisms (Enterococcus faecium, Lactobacillus acidophilus and Bifidobacterium), unique and proprietary L-form Lactobacillus, microbial sugars, enzymes and specialized proteins.

## CONVERT Users Report:

- ✔ Reduced Scours
- ✔ Increased Appetite
- ✔ Reduced Antibiotic Usage
- ✔ Reduced Mortality

### STEP 1 Day 1 - One Bolus or 15cc Gel or 5g BIG Calf Powder



#### CONVERT™ Calf Bolus

5 billion colony-forming units (CFU's)

##### Available Unit Sizes:

- 40 count jar
- 10 count foil pack



#### CONVERT™ Day-One Gel

5 billion colony-forming units (CFU's)

##### Available Unit Sizes:

- Multi-dose tube (60cc)
- Single-dose tube (15cc)



#### CONVERT™ BIG Calf Powder

5 billion colony-forming units (CFU's)

##### Available Unit Sizes:

- 2.2 lb jar - 200 dose
- 12.5 l pail - 1,100 dose

### STEP 2 Day 2 to weaning - 1g BIG Calf Powder daily

### STEP 3 When a challenge arises

During periods of challenge, increase daily dose to Challenge level for 1 to 2 feedings, evaluate calf response, then return to normal dose rate.

**CONVERT Calf Bolus Challenge Level - 1 Bolus**

**CONVERT Day-One Gel Challenge Level - 15cc**

**CONVERT BIG Calf Powder Challenge Level - 5g daily**